

University of Dundee

An isotope dilution based-targeted and non-targeted carbonyl neurosteroid/steroid profiling

Sharp, Sheila; Mitchell, Scott; Vallée, Monique; Kuzmanova, Elena; Cooper, Michelle; Belelli, Delia

Published in:
Analytical Chemistry

DOI:
[10.1021/acs.analchem.8b00055](https://doi.org/10.1021/acs.analchem.8b00055)

Publication date:
2018

Document Version
Peer reviewed version

[Link to publication in Discovery Research Portal](#)

Citation for published version (APA):

Sharp, S., Mitchell, S., Vallée, M., Kuzmanova, E., Cooper, M., Belelli, D., Lambert, J., & Huang, J. (2018). An isotope dilution based-targeted and non-targeted carbonyl neurosteroid/steroid profiling. *Analytical Chemistry*, 90(8), 5247-5255. <https://doi.org/10.1021/acs.analchem.8b00055>

General rights

Copyright and moral rights for the publications made accessible in Discovery Research Portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from Discovery Research Portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain.
- You may freely distribute the URL identifying the publication in the public portal.

Take down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Supporting Information

An isotope dilution based-targeted and non-targeted carbonyl neurosteroid/steroid profiling

Sheila Sharp¹, Scott Mitchell¹, Monique Vallée², Elena Kuzmanova¹, Michelle Cooper¹, Delia Belelli¹, Jeremy J. Lambert¹, and Jeffrey T.-J. Huang^{1*}

¹School of Medicine, University of Dundee, Dundee, DD1 9SY, UK

²Univ. Bordeaux, Neurocentre Magendie, Physiopathologie de la plasticité neuronale, INSERM U1215, F-33000, Bordeaux, France

Short title: ID-TNT-carbonyl NS profiling

Corresponding author:

*Corresponding author:

Jeffrey T.-J. Huang

Jacqui Wood Cancer Centre

School of Medicine

University of Dundee

Dundee, DD1 9SY

Tel: +44 (0)1382 386901

Email: j.t.j.huang@dundee.ac.uk

Content

Table S-1. The seed file for data processing using SIEVE.

Table S-2. The effects of extraction buffers to the measurement of five neurosteroids using ID-TNT-carbonyl NS profiling.

Table S-3. A cross method comparison between ID-TNT-NS profiling and GC-MS analysis.

Figure S-1. The effects of acute stress on cortex neurosteroids/steroids levels.

Figure S-2. The effect of finasteride on neurosteroid/steroid profiles in the cerebellum, cortex, hippocampus and hypothalamus.

Table S-1. The seed file for data processing using SIEVE.

MZ	RTStart	RTStop	Description
334.274	14.7	15	3 β ,5 α -THPROG
338.299	14.7	15	d4-3 β ,5 α -THPROG
334.274	15	15.5	3 α ,5 α -THPROG
338.299	15	15.5	d4-3 α ,5 α -THPROG
347.269	14.2	14.8	5 α -DHPROG
353.307	14.2	14.8	d6-5 α -DHPROG
345.254	14.2	14.6	Progesterone
354.310	14.2	14.6	d9-Progesterone
332.258	14.4	14.7	Pregnenolone
336.284	14.4	14.7	d4-Pregnenolone

Table S-2. The effects of extraction buffers to the measurement of five neurosteroids using ID-TNT-carbonyl NS profiling. Each condition was tested in triplicate.

The whole brains from three adult C57BL/6J mice (male, 4-6 months) were homogenized in methanol (0.5 ml / 100 mg tissue) in the presence of internal standards. The homogenate was divided into 3 parts. For each part, the composition of extraction buffer was adjusted to either 75% methanol/25% water, 1% acetic acid in methanol, or 1.8% formic acid in methanol. Each homogenate is further homogenized before subjecting to the sample preparation procedure described above. Each condition was tested in triplicate.

(pmol/g tissue)	MeOH:water 75:25 v/v	MeOH:acetic acid 99:1 v/v	MeOH:formic acid 98.2:1.8 v/v	P value (ANOVA)
3 α ,5 α -THPROG	6.5 \pm 0.7	10.5 \pm 3.5	5.0 \pm 0.0	<0.01
3 β ,5 α -THPROG	7.0 \pm 0.0	8.5 \pm 0.7	3.5 \pm 0.7	<0.01
5 α -DHPROG	33.5 \pm 0.7	22.0 \pm 0.0	21.0 \pm 4.2	<0.01
PROG	13.0 \pm 0.0	12.0 \pm 1.4	11.0 \pm 1.4	0.18
PREG	10.0 \pm 1.4	14.0 \pm 1.4	24.0 \pm 7.1	<0.01

Table S-3. A cross method comparison between ID-TNT-NS profiling and GC-MS analysis.

Two samples (one from extracts of whole brains of control mice, and the other further spiked with 5pmol/g of 3 α ,5 α -THPROG, and 3 β ,5 α -THPROG) were analyzed using ID-TNT-NS based LC-MS, or GC-MS method. Both systems were calibrated with the calibrators carefully prepared from the same source by the same operator.

Method	Neurosteroid/steroid	Sample A (pmol/g tissue)	Sample B (Sample A + 5pmol/g tissue spike-in)	% Recovery
ID-TNT- carbonyl NS profiling	3 α ,5 α -THPROG	8.0 \pm 0.3	12.1 \pm 0.3	93
	3 β ,5 α -THPROG	18.1 \pm 0.9	22.8 \pm 0.4	99
GC-MS	3 α ,5 α -THPROG	9.4 \pm 1.4	12.9 \pm 2.2	90
	3 β ,5 α -THPROG	20.3 \pm 0.8	25.5 \pm 1.6	101

Figure S-1. The effects of acute stress on cortex neurosteroids/steroids levels.

The neurosteroid/steroid profiles of cortex from adult mice (n = 4) subjected to a single exposure to raised platform for 10 min., or control animals, were analyzed using ID-TNT-carbonyl NS profiling.

- A. The box-and-whisker plots of five targeted neurosteroids/steroids in the cerebellum, cortex, hippocampus, and hypothalamus from the control (clear) and stressed mice (gray). * < 0.05 (independent t test). The plots reveal that stress produced a significant increase of $3\alpha,5\alpha$ -THPROG, $3\beta,5\alpha$ -THPROG and 5α -DHPROG in the cortex and PROG in hippocampus and hypothalamus.
- B. The base peak chromatograms of d4-THPROGs, THPROGs, corticosterone/11-deoxycortisol (m/z= 377.243), putative 5α -dihydrodeoxycorticosterone (m/z=363.264), putative tetrahydrodeoxycorticosterone (DHDOC)/hydroxy-THPROG (m/z=350.269) and putative hydroxyprogesterone/deoxycorticosterone (m/z=361.249) in the cortex samples from a control (in gray) and stressed (in red) animals. The y axis is in the same scale for both plots.
- C. Semi-quantification of corticosterone/11-deoxycortisol in four brain regions. The plot show that stress increased corticosterone/11-deoxycortisol level in the cortex. *, $p < 0.01$ (t test).
- D. A scatter plot of corticosterone/11-deoxycortisol vs putative 5α -dihydrodeoxycorticosterone levels. Each of four brain regions from each mouse are displayed. Control and stressed animals are labelled in green and red, respectively. Cere: cerebellum; Hippo: hippocampus, Hypo: hypothalamus. $Rho = 0.93$, $p < 0.001$ (Pearson correlation).

Figure S-2. The effect of finasteride on neurosteroid/steroid profiles in the cerebellum, cortex, hippocampus and hypothalamus.

Mice (n=4, each group) were treated with 50mg/kg finasteride (prepared in 20% 2-hydroxypropyl)- β -cyclodextrin in saline), 20% 2-hydroxypropyl)- β -cyclodextrin in saline, or

sham by a single s.c. injection. Two hours following the procedure mice were sacrificed and the hypothalamus, cerebellum, cortex and hippocampus dissected from each and neurosteroids were analyzed using the ID-TNT-carbonyl-NS profiling. The expression of five targeted neurosteroids/steroids in hypothalamus tissues from animals treated with finasteride, saline or sham are shown as box whisker plots. *, $p < 0.01$ vs Sham and Saline; **, $p < 0.05$ vs Saline (one-way ANOVA or Kruskal-Wallis test). Note that finasteride treatment drastically increased PREG levels in all brain regions.




Figure S-1




Figure S-2